


TAS SCAN MAPPER

Larry Kneller


1/4/07

QUICKSTART


TAS Scan Mapper displays the position in reciprocal space of the points in your scan and the location of possible spurious features in your triple axis scan. To begin open DAVE and select Planning->TAS Tools->TAS Scan Mapper



or



Two windows will appear, Scan Mapper and Scan Mapper Controls. The Scan Mapper window contains 3 plots. There is a plot of features in the scattering plane, a plot of the current scan, and a plot of the energy transfer measured. The scattering plane and the energy plot both label the scan points and all of the harmonics. The color coding is consistent between the two plots. The different harmonics can be selected for display by using the color-coded checkbox controls under the "Display Choices" tab (see below.)



The controls for the application are on a floating control base. If you delete the controls, you can restore them under the “Options” menu at the top of the main window.

The General Inputs tab contains two sub-tabs. The “Scan Definition” tab is currently implemented, however the “Setup Info” tab is implemented but its inputs are not currently used. The Setup Info tab inputs information for the resolution function, which will be implemented at a later time.

Scan Mapper controls

Scan List

General Inputs

Other Inputs

Plot Ranges

DisplayChoices

Spurion List

Scan Definition

Setup Info

Mono/Ana

D spacings

PG(002)

dm3.354

PG(002)

da3.354

Horizontal

Collimators

ALZ60

ALM60

ALA60

AL360

Vertical

Collimators

VBETZ120

VBET1120

VBET2120

VBET3120

Mosaics

ETASH10

ETASV10

ETAM30

ETAA30

SourceThermal

☐ Horizontal Focus

Instrument

Orientation

epm1

ep1

Calculate

The “Other Inputs” tab contains some important input values. Notably, the sample can material, the axis to use for the scan and energy plots, the tolerances for determining whether points in the scan are likely to give rise to spurious data, and the maximum 2-theta angle available for the spectrometer.

The screenshot shows the 'Scan Mapper controls' window with the 'Other Inputs' tab selected. The window has a title bar with standard Windows controls. Below the title bar is a tabbed interface with tabs for 'Scan List', 'General Inputs', 'Other Inputs', 'Plot Ranges', 'DisplayChoices', and 'Spurion List'. The 'Other Inputs' tab contains the following controls:


- Sample Can/Holder Material:** A text label followed by a dropdown menu showing 'Al'.
- Powder Ring:** A dropdown menu showing 'Al'.
- Qhkl values for Z-Q Plane Calculation:** Three input fields for 'Qh' (value 1), 'Qk' (value 1), and 'Ql' (value 0).
- Calculate Z-Q Plane:** A button.
- Scan Axis:** A dropdown menu showing 'k'.
- Qtol [1/Å]:** An input field with the value '0.0200000'.
- Etol [meV]:** An input field with the value '1'.
- A4max [deg]:** An input field with the value '104.000'.

The “Plot Ranges” tab gives the user control over the ranges in the three plots. The plot ranges also can be controlled directly in the plots using rubber-band box zooming.


The screenshot shows the 'Scan Mapper controls' window with the 'Plot Ranges' tab selected. The window has the same title bar and tabbed interface as the previous image. The 'Plot Ranges' tab contains three sections of controls:

- E Plot ranges:** Four input fields for 'xmin[?]' (2.90000), 'xmax[?]' (3.10000), 'dEmin[meV]' (-220.500), and 'dEmax[meV]' (220.500).
- Scattering Plane Plot ranges (1/Å):** Four input fields for 'Qx min[1/Å]' (-4.19770), 'Qx max[1/Å]' (4.19770), 'Qy min[1/Å]' (-4.19770), and 'Qy max[1/Å]' (4.19770).
- Scan Plot ranges:** Four input fields for 'xmin[?]' (2.90000), 'xmax[?]' (3.10000), 'Imin[Counts]' (0.000000), and 'Imax[Counts]' (9293.91). Below these are two checkboxes: 'x log' and 'y log', both of which are currently unchecked.


The “Display Choices” tab provides checkboxes so that the user can decide which features to display in the plots. The first set chooses Q vector display, incoherent, data display, sample can spheres, etc. The second set of checkboxes determines which types of spurious to show in the scan and energy plots. The final set of checkboxes determines which harmonics to include in the scattering plane and energy plots.




The final tab in the controls panel lists the set of spurious features along with their sources.




To define a scan you can either enter it manually via the scan definition tab in the controls or you can read in a data file. If you read in a data file, the data in the file will be displayed in the scan plot with possibly spurious features highlighted with the symbols indicating their respective sources. A text listing of the spurions will be in the spurion list tab.





All of the scans entered by the user or read in from file will appear in the Scan list in the Scan Mapper controls. Any of the scans in the list can be re-displayed by clicking on its item in the scan tree. The session can be saved by going to File->Save. This will save your work in a .spu file. Saved sessions can be recovered by selecting File->Restore Session. Note that this will restore the previous file and remove your current scans, so be sure to save your scans before restoring old work.



Finally, you can save any of the plots in either a .jpg or .ps file.

Future Features

Calculated scan of resolution-limited feature.

Resolution Ellipsoid calculations (for display and determination of spurious data)

Reciprocal Lattice calculator (Allow user input of unit cell information)

Auxiliary 3d view of reciprocal space

Output features as requested and defined by users