

QUANTUM ROTATIONS IN METHYL IODIDE: A Study Using HFBS and FANS




High Flux Backscattering Spectrometer




High Flux Backscattering Spectrometer





resolution ~ 1 μeV
energy transfers upto ±50 μeV
motions on ns timescales
momentum transfers from 0.25 Å⁻¹ to 1.75 Å⁻¹


BULK CH₃I DYNAMICS




Use 3-fold symmetric potential

CH₃I melting point = 206 K

HFBS DATA FOR CH₃I at T = 8 K




Filter Analyzer Neutron Spectrometer




measures vibrational motion
energy transfers of order 10 to 100 meV
resolution about 1 meV


FANS DATA at **T** = 12 K





HFBS DATA FOR CH₃I at T = 40 K



EISF Fitting for T = 40 K data







SUMMARY

- HFBS was used to study the rotational dynamics of CH₃I.
- Tunneling energy was determined using HFBS from which the barrier height and the transition energy for libration was calculated.
- Data from FANS was used to verify the prediction of the libration transition energy.
- Radius of the methyl group was estimated by fitting the EISF data.

ACKNOWLEDGEMENTS

•To NIST and NSF for the opportunity to experiment and study Neutron Spectroscopy and its applications

•To the NCNR Summer School Coordinators and Instructors for their assistance and guidance